

New

LIFE SCIENCE

TCI

Endoglycoceramidase

Hydrolytic Enzymes Specific to Glycosphingolipids

rEGCase I

300 mU/vial [R0240]

rEGCase II

100 mU/vial [R0242]

rEGCase I assisted by Activator II


300 mU/vial [R0241]

rEGCase II assisted by Activator II


100 mU/vial [R0243]

Endoglycoceramidase (EGCase) is a glycolipid-specific hydrolase that cleaves the glycosidic linkage between oligosaccharide and ceramide of various glycosphingolipids (GSLs).

Activity of EGCase (Ex.: Ganglioside GM1)



Hydrolysis reaction of GM1 using EGCases



Recombinant EGCase (rEGCase) is activated under a detergent-dependent condition. On the other hand, Activator II is capable of inducing activity of EGCases without any detergent reagents. By the use of Activator II, GSLs on cell surfaces of living cells could be hydrolyzed without cell disruption caused by detergent.

The products were commercialized under license from Kyushu University.

Endoglycoceramidase

Hydrolytic Enzymes Specific to Glycosphingolipids

Substrate specificities of rEGCase I and rEGCase II

| | rEGCase I | | | | rEGCase II | | | |
|-------------------|----------------|-----|-----|-----|------------|-----|-----|-----|
| Reaction time (h) | 1 | | 16 | | 1 | | 16 | |
| Enzyme (mU) | 0.4 | 4 | 0.4 | 4 | 0.4 | 4 | 0.4 | 4 |
| Substrates | Hydrolysis (%) | | | | | | | |
| Ganglio-series | | | | | | | | |
| GM1 | 89 | 100 | 100 | 100 | 82 | 100 | 93 | 100 |
| Asialo GM1 | 63 | 100 | 100 | 100 | 82 | 100 | 100 | 100 |
| Fucosyl GM1 | 54 | 92 | 92 | 96 | 20 | 38 | 27 | 52 |
| GM3 | 100 | 100 | 100 | 100 | 100 | 100 | 100 | 100 |
| GD1 | 79 | 97 | 90 | 100 | 78 | 90 | 86 | 91 |
| GT1 | 35 | 78 | 83 | 86 | 62 | 63 | 78 | 81 |
| Globo-series | | | | | | | | |
| Gb4 | 4 | 29 | 26 | 60 | 4 | 16 | 6 | 20 |
| Lacto-series | | | | | | | | |
| LacCer | 50 | 64 | 81 | 100 | 63 | 81 | 83 | 100 |
| Cerebrosides | | | | | | | | |
| GlcCer | 0 | 5 | 5 | 27 | 0 | 0 | 0 | 0 |
| GalCer | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |

Each substrate (4 nmol) was incubated with EGCase and 0.1% Triton™ X-100 at 37°C.

These data were provided by Prof. Makoto Ito and Dr. Yohei Ishibashi, Kyushu University.

References

- 1) M. Ito, T. Yamagata, *J. Biol. Chem.* **1986**, 261, 14278.
[https://doi.org/10.1016/S0021-9258\(18\)67015-2](https://doi.org/10.1016/S0021-9258(18)67015-2)
- 2) M. Ito, T. Yamagata, *J. Biol. Chem.* **1989**, 264, 9510.
[https://doi.org/10.1016/S0021-9258\(18\)60561-7](https://doi.org/10.1016/S0021-9258(18)60561-7)
- 3) Y. Ishibashi, U. Kobayashi, A. Hijikata, K. Sakaguchi, H. M. Goda, T. Tamura, N. Okino, M. Ito, *J. Lipid Res.* **2012**, 53, 2242.
<https://doi.org/10.1194/jlr.D028951>
- 4) M. Ito, Y. Ikegami, T. Tai, T. Yamagata, *Eur. J. Biochem.* **1993**, 218, 637.
<https://doi.org/10.1111/j.1432-1033.1993.tb18417.x>
- 5) M. Ito, Y. Ikegami, T. Yamagata, *Eur. J. Biochem.* **1993**, 218, 645.
<https://doi.org/10.1111/j.1432-1033.1993.tb18418.x>

Related Products

Ganglioside GM₁
TRITON™ X-100

Contact us **[G0483]**
5g / 25g **[P1775]**

For further information please refer to our website at www.TCIchemicals.com. ▶▶▶

TCI EGCase



Ordering and Customer Service

TCI AMERICA

Tel : 800-423-8616 / 503-283-1681
Fax : 888-520-1075 / 503-283-1987
E-mail : Sales-US@TCIchemicals.com

TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00
Fax : +32 (0)3 735 07 01
E-mail : Sales-EU@TCIchemicals.com

TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00
Fax : +49 (0)6196 64053-01
E-mail : Sales-DE@TCIchemicals.com

Tokyo Chemical Industry UK Ltd.

Tel : +44 (0)1865 784560
Fax : +44 (0)1865 784561
E-mail : Sales-UK@TCIchemicals.com

TCI Chemicals (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909
Fax : 044-2262 8902
E-mail : Sales-IN@TCIchemicals.com

梯希爱(上海)化成工业发展有限公司

Tel : 800-988-0390 / 021-67121386
Fax : 021-6712-1385
E-mail : Sales-CN@TCIchemicals.com

TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878
Fax : +81 (0)3-5640-8902
E-mail : globalbusiness@TCIchemicals.com

Availability, price or specification of the listed products are subject to change without prior notice. Reproduction forbidden without the prior written consent of Tokyo Chemical Industry Co., Ltd.